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\* AECOPD = acute exacerbations of chronic obstructive pulmonary disease

1. Janssens JP, Howarth-Frey C, Chevrolet JC, Abajo B, and Rochat T. Transcutaneous  $pCO_2$  to monitor noninvasive mechanical ventilation in adults - assessment of a new transcutaneous  $pCO_2$ . *Chest* 1998; 113: 768-73. 2. AARC Clinical Practice Guideline, Patient-Ventilator System Checks; Reprinted from the August 1992 issue of *Respiratory Care* *Respir Care* 1992; 37(8): 882-86. 3. Cox M, Kemp R, Anwer S *et al.* Non-invasive monitoring of  $CO_2$  levels in patients using NIV for AECOPD. *Thorax*. *bmjournals.com* 2005; Postscript; doi: 10.1136/thx.2005051664. 4. Török SS, Leuppi D *et al.* Combined oximetry-cutaneous capnography in patients assessed for long-term oxygen therapy. *Chest* 2008; 133: 1421-25.

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5. Tobias JD. Transcutaneous carbon dioxide monitoring in infants and children. *Paediatr Anaesth* 2009; 19: 434–44. 6. Johnson DC, Batool S, Dalbec B. Transcutaneous carbon dioxide monitoring in a specialized weaning unit. *Respir care* 2008 Aug; 53 (8): 1042–47. 7. Parker JC, Hernandez LA, Longenecker GL, Peevy KJ, Johnson W. Lung edema caused by high peak inspiratory pressures in dogs: role of increased microvascular filtration pressure and permeability. *Am Rev Respir Dis* 1990; 142: 321-28. 8. Dreyfuss D, Soler P, Basset G, Saumon G. High inflation pressure pulmonary edema: respective effects of high airway pressure, high tidal volume, and positive end-expiratory pressure. *Am Rev Respir Dis* 1988; 137: 1159-64. 9. Hernandez LA, Peevy KJ, Moise AA, Parker JC. Chest wall restriction limits high airway pressure-induced lung injury in young rabbits. *J Appl Physiol* 1989; 66: 2364-68. 10. Dreyfuss D, Saumon G. Ventilator-induced lung injury: lessons from experimental studies. *Am J Respir Crit Care Med* 1998; 157: 294-323. 11. American Thoracic Society, European Society of Intensive Care Medicine, Societe de Reanimation de Langue Fran-caise. International Consensus Conferences in Intensive Care Medicine: ventilator-associated lung injury in ARDS. *Am J Respir Crit Care Med* 1999; 160: 2118-24.

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12. Ahmad. Pre & post operative complications of anesthesia. The Professional 1995; 2: 1-5. 13. Gerald J. Kast. Principles & practice of point of care testing. Philadelphia: Lippincott Williams & Wilkins, 2002. 14. Hautmann H, Gamarra F et al. High Frequency Jet Ventilation in interventional fiberoptic bronchoscopy. Anesth Analg. 2000; 90: 1436-40. 15. Rowley DD, Walsh BK et al. Evaluation of a new digital transcutaneous  $tcpCO_2$  and  $SpO_2$  combination sensor and its correlation to ABG  $PaCO_2$  measurements. AARC 51st international respiratory congress, Program from the open forum no. 1. 2005; OF-05-106: 69. 16. Biro P, Layer M et al. Carbon dioxide elimination during high-frequency jet ventilation for rigid bronchoscopy. Br Journal Anaesth 2000; 80; 5: 635-37. 17. Fontana GA, Lavorini F, Boddi V et al. Comparison of dose-response slopes obtained by simultaneous assessment of changes in FEV1 and transcutaneous oxygen partial pressure during methacholine challenges in asthmatic patients. J Asthma 1995; 32, 4: 301-407. 18. Fontana GA, Cardellicchio S, Camiciottoli G et al. Changes in transcutaneous oxygen partial pressure as an index of response to inhaled methacholine in asthmatic patients. Chest 1993; 103,5: 1375-80. 19. Palmisano BW, Severinghaus JW. Transcutaneous  $pCO_2$  and  $pO_2$  A multicenter study of accuracy. J Clin Monit 1990; 6: 189-95. 20. Bradley S, Solin P, Wilson J et al. Hypoxemia and hypercapnia during exercise and sleep in patients with cystic fibrosis. Chest 1999; 116(3): 647-54.

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